

Vaistų koncentracijos ribinis diapazonas	OPI/MOP40		K225		MTD30		OXY20		MDMA50		BZO50		BAR50	
	-	+	-	+	-	+	-	+	-	+	-	+	-	+
0% Riba	30	0	30	0	30	0	30	0	30	0	30	0	30	0
-50% Riba	30	0	30	0	30	0	30	0	30	0	30	0	30	0
-25% Riba	27	3	26	4	25	5	27	3	26	4	25	5	25	5
Riba	13	17	15	15	15	15	20	10	19	11	13	17	18	12
+25% Riba	7	23	3	27	7	23	4	26	6	24	4	26	8	22
+50% Riba	0	30	0	30	0	30	0	30	0	30	0	30	0	30
+300% Riba	0	30	0	30	0	30	0	30	0	30	0	30	0	30

Vaistų koncentracijos ribinis diapazonas	THC20		THC50		OXY50		LSD10		BZO30		BZO10		BUP10		PPX50	
	-	+	-	+	-	+	-	+	-	+	-	+	-	+	-	+
0% Riba	30	0	30	0	30	0	30	0	30	0	30	0	30	0	30	0
-50% Riba	30	0	30	0	30	0	30	0	30	0	30	0	30	0	30	0
-25% Riba	27	3	27	3	27	3	26	4	25	5	26	4	26	4	26	4
Riba	12	18	14	16	20	10	16	14	13	17	19	11	14	16	16	14
+25% Riba	8	22	9	21	4	26	7	23	4	26	6	24	10	20	6	24
+50% Riba	0	30	0	30	0	30	0	30	0	30	0	30	0	30	0	30
+300% Riba	0	30	0	30	0	30	0	30	0	30	0	30	0	30	0	30

Vaistų koncentracijos ribinis diapazonas	MQL300		THC40		CAR300		EDDP50		ABP10		UR-14425		ZOP50		BUP5		FYL10	
	-	+	-	+	-	+	-	+	-	+	-	+	-	+	-	+	-	+
0% Riba	30	0	30	0	30	0	30	0	30	0	30	0	30	0	30	0	30	0
-50% Riba	30	0	30	0	30	0	30	0	30	0	30	0	30	0	30	0	30	0
-25% Riba	27	3	27	3	28	2	27	3	25	5	28	2	27	3	26	4	26	4
Riba	15	15	14	16	16	14	14	16	15	15	15	17	13	14	16	15	15	15
+25% Riba	4	26	9	21	3	27	4	26	4	26	3	27	4	26	10	20	3	27
+50% Riba	0	30	0	30	0	30	0	30	0	30	0	30	0	30	0	30	0	30
+300% Riba	0	30	0	30	0	30	0	30	0	30	0	30	0	30	0	30	0	30

Analitinis specifškumas
Toliau pateiktoje lentelėje išvardytos junginių koncentracijos (ng/mL), viršijančios Multi-Drug greitojo testo taurelės skirtos AMP/MET/COC/OPI/MOP/THC/PCP/MTD/MDMA/BZO/OXY/COT/K2/KET/BAR/BUP/6-MAM/TML/FYL/CFYL/MDPV/ α -PVP/LSD/PPX/MQL/CAR/EDDP/ABP/UR-144/ZO P nustatę teigiamus rezultatus per 3–10 minučių skaitymo laiką.

Compound	ng/mL	Compound	ng/mL
Amphetamine (AMP 50)			
d-Amphetamine	50	p-Hydroxyamphetamine	100
d/l-Amphetamine	100	(+)-3,4-Methylenedioxyamphetamine (MDA)	100
β -Phenylethylamine	25,000	l-Amphetamine	25,000
Tryptamine	12,500	Methoxyphenamine	12,500
Methamphetamine (MET 50)			
d-Methamphetamine	50	l-Phenylephrine (R)-(-)-Phenylephrine	6,250
Fenfluramine	60,000	Procaine	2,000
p-Hydroxymethamphetamine	400	(1R,2S)-(-)-Ephedrine	400
Methoxyphenamine	25,000	Ephedrine	400
Mephentermine	1,500	Benzphetamine	25,000
3,4-Methylenedioxyamphetamine (MDMA)	50		
Marijuana (THC 12)			
11-nor- Δ^9 -THC-9 COOH	12	Δ^9 -THC	10,000
Cannabinol	12,500	11-nor- Δ^8 -THC-9 COOH	12
Δ^8 -THC	6,000		
Marijuana (THC 50)			
11-nor- Δ^9 -THC-9 COOH	50	Δ^9 -THC	40,000
Cannabinol	50,000	11-nor- Δ^8 -THC-9 COOH	40
Δ^8 -THC	25,000		
Marijuana (THC 20)			
11-nor- Δ^9 -THC-9 COOH	20	Δ^9 -THC	17,000
Cannabinol	20,000	11-nor- Δ^8 -THC-9 COOH	15
Δ^8 -THC	10,000		
Marijuana (THC 40)			
11-nor- Δ^9 -THC-9 COOH	40	Δ^9 -THC	35,000
Cannabinol	40,000	11-nor- Δ^8 -THC-9 COOH	35
Δ^8 -THC	20,000		

Cocaine (COC 20)			
Benzoylecgonine	20	Egonine	1,500
Cocaine	20	Egonine methyl ester	12,500
Cocacethylene	30		
Opiates (OPI/MOP 40)			
Morphine	40	Norcodeine	6,250
Codeine	25	Normorphine	25,000
Ethylmorphine	25	Nalorphine	10,000
Hydromorphone	100	Oxymorphone	25,000
Hydrocodone	100	Thebaine	2,000
Levorphanol	400	Diacetylmorphine (Heroin)	50
Oxycodone	25,000	6-Monoacetylmorphine	25
Morphine 3- β -D-Glucuronide	50		
Phencyclidine (PCP 10)			
Phencyclidine	10	4-Hydroxyphencyclidine	2,500
Oxycodone (OXY 20)			
Oxycodone	20	Hydromorphone	10,000
Oxymorphone	40	Naloxone	5,000
Levorphanol	10,000	Naltrexone	5,000
Hydrocodone	1,500		
Oxycodone (OXY 50)			
Oxycodone	50	Hydromorphone	20,000
Oxymorphone	100	Naloxone	12,500
Levorphanol	25,000	Naltrexone	12,500
Hydrocodone	3,750		
Cotinine (COT 20)			
(-)-Cotinine	20	(-)-Nicotine	300
Synthetic Marijuana (K2-25)			
JWH-018 5-Pentanoic acid metabolite	25	JWH-018 5-Hydroxypentyl metabolite	250
JWH-073 4-butanoic acid metabolite	25	JWH-073 4-Hydroxybutyl metabolite	250
JWH-018 4-Hydroxypentyl metabolite	200		
Benzodiazepines (BZO 50)			
Alprazolam	25	Flunitrazepam	25
a-hydroxyalprazolam	250	(\pm) Lorazepam	500
Bromazepam	130	RS-Lorazepamglucuronide	25
Chlordiazepoxide	130	Midazolam	1,000
Clobazam	25	Nitrazepam	25
Clonazepam	65	Norchlordiazepoxide	25
Clorazepatedipotass	65	Nordiazepam	130
Delorazepam	130	Oxazepam	50
Desalkylflurazepam	25	Temazepam	25
Diazepam	250	Triazolam	500
Estazolam	1,000		
Benzodiazepines (BZO 30)			
Alprazolam	15	Flunitrazepam	15
a-hydroxyalprazolam	150	(\pm) Lorazepam	300
Bromazepam	75	RS-Lorazepamglucuronide	15
Chlordiazepoxide	75	Midazolam	600
Clobazam	15	Nitrazepam	15
Clonazepam	40	Norchlordiazepoxide	15
Clorazepatedipotass	40	Nordiazepam	75
Delorazepam	75	Oxazepam	30
Desalkylflurazepam	15	Temazepam	15
Diazepam	150	Triazolam	300
Estazolam	600		
Benzodiazepines (BZO 20)			
Alprazolam	10	Flunitrazepam	10
a-hydroxyalprazolam	100	(\pm) Lorazepam	200
Bromazepam	50	RS-Lorazepamglucuronide	10
Chlordiazepoxide	50	Midazolam	400
Clobazam	10	Nitrazepam	10
Clonazepam	25	Norchlordiazepoxide	10
Clorazepatedipotass	25	Nordiazepam	50
Delorazepam	50	Oxazepam	20
Desalkylflurazepam	10	Temazepam	10
Diazepam	100	Triazolam	200
Estazolam	400		
Benzodiazepines (BZO 10)			
Alprazolam	10	Flunitrazepam	10
a-hydroxyalprazolam	80	(\pm) Lorazepam	150
Bromazepam	40	RS-Lorazepamglucuronide	10
Chlordiazepoxide	40	Midazolam	300
Clobazam	10	Nitrazepam	10
Clonazepam	20	Norchlordiazepoxide	10
Clorazepatedipotass	20	Nordiazepam	40
Delorazepam	40	Oxazepam	10
Desalkylflurazepam	10	Temazepam	10
Diazepam	80	Triazolam	150
Estazolam	300		
Methadone (MTD 30)			
Methadone	30	LAAM	200
Disopyramide	400	Doxylamine	12,500
(+)-Chlorpheniramine	6,250	Nor-LAAM	12,500
Methylenedioxyamphetamine (MDMA 50)			
(\pm) 3,4-Methylenedioxyamphetamine HCl (MDMA)	50	3,4-Methylenedioxyethylamphetamine (MDE)	30
(\pm) 3,4-Methylenedioxyamphetamine HCl (MDA)	300	l-Methamphetamine	25,000
Ketamine (KET 50)			
Ketamine	50	Mephentermine	1,250
Tetrahydrozoline	20	Phencyclidine	625
Benzphetamine	1,250	(1R,2S)-(-)-Ephedrine	5,000
d-Methamphetamine	1,250	Promazine	1,250
(+)-Chlorpheniramine	1,250	4-Hydroxyphencyclidine	2,500
l-Methamphetamine	2,500	Promethazine	1,250
Clomidine	5,000	Levorphanol	2,500
Methoxyphenamine	625	Thioridazine	2,500
Disopyramide	625	MDE	2,500
d-Norpropoxyphene	625	Meperidine	1,250
EDDP	2,500	Dextromethorphan	75
Pentazocine	1,250	(+)-3,4-Methylenedioxyamphetamine (MDMA)	5,000
Barbiturates (BAR 50)			
Amobarbital	800	Alphenol	100
5,5-Diphenylhydantoin	1,500	Aprobarbital	80
Allobarbitol	100	Butabarbitol	40
Barbital	1,500	Butalbitol	1,500
Talbutal	40	Butethal	90
Cyclopentobarbital	5,000	Phenobarbital	50
Pentobarbital	1,500	Secobarbital	50

Buprenorphine (BUP 10)			
Buprenorphine	10	Norbuprenorphine	50
Buprenorphine 3-D-Glucuronide	50	Norbuprenorphine 3-D-Glucuronide	100
Buprenorphine (BUP 5)			
Buprenorphine	5	Norbuprenorphine	25
Buprenorphine 3-D-Glucuronide	25	Norbuprenorphine 3-D-Glucuronide	50
Tramadol (TML 30)			
n-Desmethyl-cis-tramadol	60	o-Desmethyl-cis-tramadol	3,000
Cis-tramadol	30	Phencyclidine	30,000
Procyclidine	30	d,l-O-Desmethylvenlafaxine	15,000
6-Monoacetylmorphine (6-MAM 10)			
6-Monoacetylmorphine	10	Morphine	100,000
Fentanyl (FYL 50)			
Alfentanyl	1,500,000	Buspirone	37,500
Fenfluramine	125,000	Fentanyl	50
Norfentanyl	10	Sufentanyl	125,000
Fentanyl (FYL 20)			
Alfentanyl	600,000	Buspirone	37,500
Fenfluramine	50,000	Fentanyl	20
Norfentanyl	8	Sufentanyl	50,000
Fentanyl (FYL 10)			
Alfentanyl	300,000	Buspirone	20,000
Fenfluramine	25,000	Fentanyl	10
Norfentanyl	8	Sufentanyl	25,000
Carfentanyl (CFYL 50)			
Carfentanyl	50	Fentanyl	25
Sufentanil	300	(±)cis-3-Methylfentanyl	50,000
Ramifentanil	500	Butylfentanyl	200
3,4-methylenedioxypropylvalerone (MDPV 300)			
3,4-methylenedioxypropylvalerone	300		
alpha-Pyrrolidinovaleferone (α-PVP 300)			
alpha-Pyrrolidinovaleferone	300		
Lysergic Acid Diethylamide (LSD 10)			
Lysergic Acid Diethylamide	10		
Propoxyphene (PPX 50)			
d-Norpropoxyphene	50	d-Propoxyphene	50
Methaqualone (MQL 300)			
Methaqualone	300		
Carisoprodol (CAR 300)			
Carisoprodol	300		
2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP 50)			
2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP)	50		
AB-PINACA (ABP 10)			
AB-PINACA	10	UR-144 4-hydroxypentyl	10,000
AB-PINACA 5-Pentanoic	10	APINACA 5-hydroxypentyl	10,000
AB-PINACA 5-hydroxypentyl	10	AB-FUBINACA	10
ADB-PINACA N-(5-hydroxypentyl)	30	ADB-PINACA Pentanoic Acid	10
AB-PINACA 4-hydroxypentyl	10,000	5-fluoro AB-PINACAN-(4-hydroxypentyl)	30
UR-144 5-hydroxypentyl	10,000	5-fluoro AB-PINACA	25
UR-144 5-Pentanoic		AB-CHMINACA	100
UR-144			
UR-144 5-Pentanoic acid	25	UR-144 4-hydroxypentyl	10,000
5-fluoro AB-Pinaca N-(4- hydroxypentyl)	10,000	ADB-PINAC N-(4- hydroxypentyl)	>10,000
UR-144 5-hydroxypentyl	5,000	AB-PINACA 4-hydroxypentyl	>10,000
XLR-11 4-hydroxypentyl	2,000		
Zopiclone (ZOP 50)			
Zopiclone-N-oxide	50	Zopiclone	50

Kryžminis reaktyvumas

Buvo atliktas tyrimas, siekiant nustatyti bandomo kryžminį reaktyvumą su junginiais, įpiltais į PBS žaliava be narkotikų. Šie junginiai neparodė klaidingai teigiamų kelių narkotikų greitojo testo taurelės rezultatų, kai buvo tirama koncentracija iki 100 µg/mL.

Acetaminophen	d/l-Chlorpheniramine	Sulfamethazine
N-Acetylprocainamide	Chloroquine	Tetracycline
Aminopyrine	Clonidine	Tetrahydrocortisone 3 (β-D-glucuronide)
Ampicillin	l-Cotinine	Thioridazine
Apomorphine	Deoxycorticosterone	Tolbutamide
Atropine	Diclofenac	Trifluoperazine
Benzoic acid	Digoxin	d/l-Tryptophan
d/l-Brompheniramine	l-ψ-Ephedrine	Uric acid
Chloral-hydrate	Estrone-3-sulfate	Ketoprofen
Chlorothiazide	l(-)-Epinephrine	Loperamide
Chlorpromazine	Fenoprofen	Meprobamate
Cholesterol	Gentisic acid	Nalidixic acid
Cortisone	Hydralazine	Niacinamide
Creatinine	Hydrocortisone	Norethindrone
Dextromethorphan	p-Hydroxytyramine	Noscapine
Diflunisal	Iproniazid	Oxalic acid
Diphenhydramine	Isoxsuprine	Oxymetazoline
β-Estradiol	Labetalol	Penicillin-G
Ethyl-p-aminobenzoate	Meperidine	Perphenazine
Erythromycin	Methylphenidate	Trans-2-phenylcyclopropylaminehydrochloride
Furosemide	Naproxen	Prednisolone
Hemoglobin	Nifedipine	d/l-Propranolol
Hydrochlorothiazide	d/l-Octopamine	d-Pseudoephedrine
o-Hydroxyhippuric acid	Oxolinic acid	Quinine
Ibuprofen	Papaverine	Ranitidine
d/l-Isoproterenol	Pentazocine hydrochloride	Serotonin
Acetophenetidin	Phenelzine	Sulindac
Acetylsalicylic acid	Phenylpropanolamine	Tetrahydrocortisone 3-acetate
Amoxicillin	Prednisone	Thiamine
l-Ascorbic acid	Quinacrine	d/l-Tyrosine
Aspartame	Quindine	Triamterene
Benzoic acid	Salicylic acid	Trimethoprim
Benzphetamine	Zomepirac	Tyramine
Caffeine	Chloramphenicol	Verapamil

ŠALTINIAI

- Moolchan, E., et al. "Saliva and Plasma Testing for Drugs of Abuse: Comparison of the Disposition and Pharmacological Effects of Cocaine", Addiction Research Center, IRP, NIDA, NIH, Baltimore, MD. As presented at the SOFT-TIAFT meeting October 1998.
- Kim, I. et al. "Plasma and oral fluid pharmacokinetics and pharmacodynamics after oral codeine administration", *ClinChem*, 2002 Sept.; 48 (9), pp 1486-96.
- Schramm, W. et al. "Drugs of Abuse in Saliva: A Review," *J Anal Tox*, 1992 Jan-Feb; 16 (1), pp 1-9
- McCarron, MM, et al. "Detection of Phencyclidine Usage by Radioimmunoassay of Saliva," *J Anal Tox*, 1984 Sep-Oct.; 8 (5), pp 197-201.
- International Narcotics Control Board. Report of the International Narcotics Control Board for 2009[R]. New York: UN, 2010
- Lane JC, Tension MB, Lawless ST, et al. Movement disorder after withdrawal of fentanyl infusion. *J Pediatr*, 1991, 119 (4): 649-651
- Dominguez KD, Lomako DM, Katz RW, et al. Opioid withdraw in critically ill neonates. *Ann Pharmacotherm*, 2003, 37 (4): 473-477
- European Monitoring Centre for Drugs and Drug Addiction. Annual Report 2009[R]. Lisbon: EMCDDA, 2010
- Stanley, Theodore H.; Egan, Talmage D.; Aken, Hugo Van (February 2008). "A Tribute to Dr. Paul A. J. Janssen: Entrepreneur Extraordinaire, Innovative Scientist, and Significant Contributor to Anesthesiology". *Anesthesia & Analgesia*. 106 (2): 451-462.
- Vos, V. De (22 July 1978). "Immobilisation of free-ranging wild animals using a new drug". *Veterinary Record*. 103 (4): 64-68. Dominguez KD, Lomako DM, Katz RW, et al. Opioid withdraw in critically ill neonates. *Ann Pharmacotherm*, 2003, 37 (4): 473-477
- Mouteney, Jane; Giraudon, Isabelle; Denissov, Gleb; Griffiths, Paul (July 2015). "Fentanyls: Are we missing the signs? Highly potent and on the rise in Europe". *International Journal of Drug Policy*. 26 (7): 626-631.
- Glass, IB. *The International Handbook of Addiction Behavior*. Routledge Publishing, New York, NY. 1991, 216
- Baselt RC. *Disposition of Toxic Drugs and Chemicals in Man*. 6th Ed.
- Waugh; et al. (2013). "Deaths Involving the Recreational Use of α-PVP (α-pyrrolidinopentiophenone)" (PDF). AAFS Proceedings. Abstract K16.
- Glass, IB. *The International Handbook of Addiction Behavior*. Routledge Publishing, New York, NY. 1991, 216
- Kamin I, Shaskan D (1959). "Death due to massive overdose of meprobamate". *Am J Psychiatry*. 115 (12): 1123-1124. doi:10.1176/ajp.115.12.1123-a. PMID 13649976.
- Hollister LE (1983). "The pre-benzodiazepine era". *J Psychoactive Drugs*. 15 (1-2): 9-13. doi:10.1080/02791072.1983.10472117. PMID 6350551.
- Gaillard Y, Billault F, Pepin G (1997). "Meprobamate overdose: a continuing problem. Sensitive GC-MS quantitation after solid phase extraction in 19 fatal cases". *Forensic Sci. Int*. 86 (3): 173-180. doi:10.1016/S0379-0738(97)02128-2. PMID 9180026.
- Allen MD, Greenblatt DJ, Noel BJ (1977). "Meprobamate overdose: a continuing problem". *Clin Toxicol*. 11 (5): 501-515. doi:10.3109/15563657708988216. PMID 608316.
- Kintz P, Tracqui A, Mangin P, Lugnier AA (1988). "Fatal meprobamate self-poisoning". *Am J Forensic Med Pathol*. 9 (2): 139-140. doi:10.1097/0000433-198806000-00009. PMID 3381792.
- Eeckhout E, Huyghens L, Loeff B, Maes V, Sennesael J (1988). "Meprobamate poisoning, hypotension and the Swan-Ganz catheter". *Intensive Care Med*. 14 (4): 437-438. doi:10.1007/BF00262904. PMID 3403779.
- Lhoste F, Lemaire F, Rapin M (1977). "Treatment of hypotension in meprobamate poisoning". *N Engl J Med*. 296 (17): 1004. doi:10.1056/NEJM197704282961717. PMID 846530.
- Bedson H (1959). "Coma due to meprobamate intoxication. Report of a case confirmed by chemical analysis". *Lancet*. 273 (1): 288-290. doi:10.1016/S0140-6736(59)90209-0. PMID 13632000.
- Uchiyama, N.; Matsuda, S.; Wakana, D.; Kikura-Hanajiri, R.; Goda, Y. (2012). "New cannabimimetic indazole derivatives, N-(1-amino-3-methyl-1-oxobutan-2-yl)-1-pentyl-1H-indazole-3-carboxamide (AB-PINACA) and N-(1-amino-3-methyl-1-oxobutan-2-yl)-1-(4-fluorobenzyl)-1H-indazole-3-carboxamide (AB-FUBINACA) identified as designer drugs in illegal products". *Forensic Toxicology*. 31: 93-100. doi:10.1007/s11419-012-0171-4.
- AB-PINACA". *Cayman Chemical*. Retrieved 25 June 2015.
- Banister, Samuel D.; Moir, Michael; Stuart, Jordyn; Kevin, Richard C.; Wood, Katie E.; Longworth, Mitchell; Wilkinson, Shane M.; Beinart, Corinne; Buchanan, Alexandra S.; Glass, Michelle; Connor, Mark; McGregor, Iain S.; Kassiou, Michael (2015). "Pharmacology of Indole and Indazole Synthetic Cannabinoid. Designer Drugs AB-FUBINACA, ADB-FUBINACA, AB-PINACA, ADB-PINACA, 5F-AB-PINACA, 5F-ADB-PINACA, ADBICA, and 5F-ADBICA". *ACS Chemical Neuroscience*. 6 (9): 1546-59. doi:10.1021/acs.cchem.5b00100.
- Volpicellim, Joseph R., M.D., Ph.D.: *Alcohol Dependence: Diagnosis, Clinical Aspects and Psychosocial Causes*. Substance Abuse Library, University of Pennsylvania, 1997.
- Jones, A.W.: *Inter- and intra individual variations in the saliva/blood alcohol ratio during ethanol metabolism in man*. *Clin. Chem*. 25, 1394-1398, 1979.
- Baselt RC. *Disposition of Toxic Drugs and Chemicals in Man*. 2nd Ed. Biomedical Publ., Davis, CA. 1982; 488

Simboliai

	Peržiūrėkite naudojimo instrukcijas		Yra pakankamai <n> testams		Temperatūros riba
	In vitro diagnostinis medicinos prietaisas		Partijos kodas		Katalogo numeris
	Igalios atstovas Europos bendrijoje / Europos Sąjungoje		Tinkamumo laikas		Nenaudoti pakartotinai
	Nenaudokite, jei pakuotė pažeista ir vadovaukitės naudojimo instrukcijomis		Gamintojas		Dėmesio

ACRO BIOTECH, Inc.
 4650 Arrow Highway,
 Suite D-6 Montclair, CA 91763, U.S.A.
 Tel: +1 (909) 541-5085 www.acrobiotech.com



MedNet EC-REP GmbH
 Borkstrasse 10
 48163 Muenster
 Germany

Lietuvos platintojas
 MB Euromedika
 T. Ševčenkos 16, korp.2
 LT-03111 Vilnius, Lietuva
 Tel. +370 5 2151418
www.diagnostikostestai.lt

